

# MULTICUT EASY & EASY PLUS

## Metallographical Sample Cutting Machine

Multicut Easy is a manual cutting machine used in metallographical sample preparation which allows you to obtain a high quality sample surface. The Multicut Easy Plus option, on the other hand, allows you to precisely cut your large samples in one go, thanks to its movable table.

- | Powerful motor specially produced for cutting operations
- | Special anti-corrosion coating for inner cutting chamber
- | Fast cooling and easy cleaning thanks to the powerful water pump
- | IP67 certified control panel that is not affected by water and dust



## / SPECIFICATIONS

**Ø305** mm

Cut-off Wheel

**Ø115** mm

Max. Cutting  
Capacity

**3** kW

Motor Power

## / PRODUCT DETAILS



## / TECHNICAL SPECIFICATIONS

Motor Rotation Speed

**2800** rpm

Table Dimensions

**296x286** mm

T-Slot Dimension

**10** mm

Disc Kit

**32 vs 12.7** mm

Water Pump

**30** Liter/min

Table Movement Range

**226** mm\*

Dimensions (w x d x h)

**676x696x660** mm

Weight

**126** kg

Power Supply

**380V**  
**60/60** Hz

\*Multicut Easy Plus