

MULTICUT SENSE

Automatic Precision Cut-off Machine

Multicut Sense is an automatic cutting device used in metallographic sample preparation and you can obtain the most sensitive surfaces. With its 4.3" color touch screen, all parameters are easy to identify and control.

- | Precise table movement with servo motor system
- | Powerful motor specially produced for cutting operations
- | Special anti-corrosion coating
- | Fast cooling and easy cleaning thanks to the powerful water pump



/ SPECIFICATIONS

Ø200 mm

Cut-off Wheel

Ø74 mm

Max. Cutting
Capacity

3 kW

Motor Power

/ PRODUCT DETAILS



/ TECHNICAL SPECIFICATIONS

Motor Rotation Speed

600-2800 rpm

Table Dimensions

296x286 mm

Program Memory

16 Records

Disc Kit

12.7/22 mm

Cutting Modes

Straight Cut
Balanced Cut

Variable Table Speed

0.6 mm/s

Dimensions (w x d x h)

670x696x460 mm

Weight

126 kg

Power Supply

380V
50/60 Hz