

MULTIMOUNT EASY

Hot Mounting Machine

The Multimount Easy is an automatic mounting device for metallographic sample preparation. It eliminates some routine maintenance costs with its pneumatic pressure system. Meet the fastest molding times in its class thanks to its special mold design.

- | Fast results with powerful heating and cooling system
- | Quarter turn lockable safety cover system
- | Method memory for different bakelite powders
- | Ease of use: Touch panel control over icons



/ SPECIFICATIONS

220 °C

Maximum
temperature

1950 daN

Pressure

25 Records

Program Memory

/ PRODUCT DETAILS



/ TECHNICAL SPECIFICATIONS

Mount Options

Ø25-30-40mm

Mounting Durations

**Ø25 ~ 4 min
Ø30 ~ 5 min
Ø40 ~ 8 min**

Control

Color Touch Panel Screen

Pressure System

Pneumatic

Cooling Options

**Fast
Medium
Slow**

Temperature Sensor

PT 100

Dimensions (w x d x h)

370x646x430mm

Weight

35kg

Power Supply

**200-240V
50/60Hz**