

MULTIPOL EASY & EASY PLUS

Manual Grinder - Polisher

Multipol Easy is a manual grinding and polishing device that can offer you the best quality in the most economical way for your routine metallographic sample preparation processes.

Multipol Easy Plus, on the other hand, is a manual grinding and polishing device with a touch control panel that will allow you to obtain perfect surfaces even in inexperienced hands thanks to its pressure control display.

- | Single or dual disc design
- | Easy cleaning thanks to extendable and flow-adjustable taps^^
- | Quick setting for 150 and 300 rpm
- | Excellent results even in manual preparation with the display of applied pressure in Newtons^^



^^Multipol Easy Plus

/ SPECIFICATIONS

50-650 RPM

Rotation Speed

AC/0.75 KW

Motor (disc)

IP67

Control Panel

/ PRODUCT DETAILS



/ TECHNICAL SPECIFICATIONS

Working Disc Options

Ø200/260/300mm

Disc Rotation Direction

**Single Direction
Both Directions**

Control System

**Analog[^]
Color touch panel^{^^}**

Dimensions (w x d x h)

585x695x260mm *
760x695x260mm **

Weight

46kg *
76kg **

Power Supply

200-240V
50/60Hz

*Single Disc

**Double Disc

[^]Multipol Easy

^{^^}Multipol Easy Plus