

# MULTIPOL SENSE

## Fully Automatic Grinder-Polisher

Multipol Sense is a fully automatic grinding and polishing device that offers the most innovative solutions in all processes of metallographic sample preparation, thanks to its superior technological equipment.

- | Single or dual disc design
- | Easy cleaning thanks to extendable and flow-adjustable taps
- | Automatic dosing system
- | Program memory for sample preparation methods of different materials
- | Self cleaning basin system



## / SPECIFICATIONS

**50-650** RPM

Disc Options

**5** Pumps

Automatic Dosing

**0-200** RPM

Head Speed

**15** Record

Program Memory


**10-60**<sup>N</sup>

Pressure Force

**28** ml/min

Max. Dosing Rate

## / PRODUCT DETAILS



## / TECHNICAL SPECIFICATIONS

Disc Options

**Ø200/250/300**mm

Sample Capacity

**6xØ40-25**mm  
**3xØ50**mm

Control

**Color touch panel**

Rotation direction

**Both direction**

Motor (disk)

**AC/0.75**kW

Motor (head)

**Servo/0.75**kW

Dimensions ( w x d x h)

**585x695x720**mm \*  
**760x695x720**mm \*\*

Weight

**95**kg \*  
**120**kg \*\*

Power Supply

**200-240**V  
**50/60**Hz

\*Single Disc ————— \*\*Double Disc